

## EXPERIMENT 15A

---

### QUALITATIVE ANALYSIS: GROUP I ( $\text{Ag}^+$ , $\text{Hg}_2^{2+}$ , $\text{Pb}^{2+}$ )

Carry out the procedures on a “known” solution that contains all of the Group I ions. Then obtain an **unknown** from your instructor and use the same procedures to determine the Group I ion(s) present.

#### PROCEDURES:

1. Add 10 drops of 6.0 M HCl to 2.0 mL of the sample. A white precipitate will form if Group I ions are present.




Centrifuge the mixture and add another drop of HCl to verify complete precipitation. If precipitation is complete proceed to the next step. If not, add 2 more drops of 6.0 M HCl, stir, centrifuge, and again test for completeness of precipitation. The supernatant must be free of Group I ions before moving on to the next step.

2. Decant and discard the supernatant. Add 2 mL distilled water to the precipitate and heat for 3 minutes in a hot water bath. Stir occasionally while heating. Any lead present in the precipitate (in the form of  $\text{PbCl}_2$ ) will dissolve in the water as it becomes hot.



Toward the end of the heating period allow the remaining precipitate to settle to the bottom of the test tube giving a clear solution above. Turn off the heat if necessary to avoid agitation of the precipitate. When the supernatant is completely clear, draw off as much as possible into a pipette without picking up any of the precipitate and place the supernatant (still hot) in a separate test tube.

3. Allow the hot supernatant from Step 2 to cool, then add 2-3 drops of 1 M  $\text{K}_2\text{CrO}_4$ . A yellow precipitate confirms the presence of lead ion. A yellow precipitate means lead is present.



**Caution:** Some white precipitate of  $\text{PbCl}_2$  once more, could form as the solution cools. Look for yellow precipitate.


4. To the remaining precipitate of Step 2, add 10 drops of 4 M aqueous ammonia,  $\text{NH}_3\text{(aq)}$ . A black or dark grey residue (precipitate) confirms the presence of mercury.



5. Centrifuge the mixture from Step 4 and decant the clear supernatant into a separate test tube. Add 4 M  $\text{HNO}_3$  until the solution is acidic (pH = 5 or less). The formation of a white precipitate (cloudiness) confirms the presence of silver.



## EXPERIMENT 15A–GROUP I ANALYSIS–FLOW CHART



NAME \_\_\_\_\_ SECTION \_\_\_\_\_

EXPERIMENT 15A – GROUP I  
REPORT SHEET

UNKNOWN NUMBER: \_\_\_\_\_ ANALYSIS \_\_\_\_\_

**YOUR FLOW CHART: (Make the chart complete by showing all steps, even washings, all reagents at proper concentrations and their amounts, color of precipitates, and any other information that was interesting and/or important.)**